

UK Biobank

Imaging modality DXA

Version 1.0

<http://www.ukbiobank.ac.uk/>

23rd Sept 2015



This document details the procedure for the DXA performed at an Imaging assessment centre for UK Biobank.

Contents

| | |
|-------------------------------------|---|
| 1. Introduction..... | 1 |
| 2. Equipment Preparation | 2 |
| 3. Participant Preparation..... | 2 |
| 4. DXA modality measurement..... | 2 |
| 5. Download/transfer of Data? | 3 |
| 6. QC protocol..... | 3 |
| 7. Data collected | 3 |

1. Introduction

1.1: This manual details the procedure for the DXA measurement at a UK Biobank Imaging Assessment Centre.

Table 1: Sequence of assessment visit

| | Visit station | Assessments undertaken |
|---|----------------------------|---|
| 1 | Reception | <ul style="list-style-type: none"> Welcome & registration Generating a USB key for Participants |
| 2 | Eligibility Section | <ul style="list-style-type: none"> Eligibility questionnaire Consent |
| 3 | Imaging scans | <ul style="list-style-type: none"> Cardiac MRI scan Body MRI scan Brain MRI scan DXA scan of whole body, bone and joint 12-lead ECG Carotid ultrasound |
| 4 | Touchscreen | <ul style="list-style-type: none"> Touchscreen questionnaire Hearing Test Cognitive function tests |
| 5 | Interview & blood pressure | <ul style="list-style-type: none"> Interviewer questionnaire Blood pressure measurement Measurement of arterial stiffness |
| 6 | Physical measurements | <ul style="list-style-type: none"> Height (Standing and Sitting) Hip & Waist measurement Weight and Bio-impedance measurement Hand-grip strength Ultrasound Bone Densitometry Spirometry (Lung function test) |
| 7 | Sample collection & exit | <ul style="list-style-type: none"> Blood, urine and saliva sample collection |

1.2: At the start of their visit, each participant is issued with a USB key at the Reception station. This contains Participant ID, name, date of birth and gender. As the participant progresses between stations the USB key acts as an identifying token. The USB key is encrypted so can only be read by assessment centre computers. None of the participant's

test data is transferred to the USB key. At the end of the assessment visit all identifying data on the USB key is removed.

1.3: This procedure is performed by a radiographer who has received suitable training and has been granted the relevant module permissions.

2. Equipment Preparation

The iDXA instrument is calibrated to a manufacturer's phantom (GE-Lunar, Madison, WI) and undergoes a daily QC procedure. Periodic assessment of QC reports is undertaken to ensure long-term stability.

3. Participant Preparation

The participant is prepared as per point 1 above, and the DXA procedure is explained. An explanation of the physical positions required and the nature of each individual scan is explained as the assessment proceeds. Scan acquisition may be repeated as appropriate, as detailed below.

4. DXA modality measurement (GE-Lunar iDXA)

- *Whole body:* The participant is asked to lie flat on their back on the DXA couch for the whole body scan.
- *Lumbar spine:* The participant's lower legs are placed on a polystyrene block, bringing the hips and knees to 80 degrees flexion. The lumbar spine is then scanned. Up to three attempts at correct placement of the starting position are permitted.
- *Hip:* The participant's foot, with leg straight, is strapped against a support to ensure correct orientation of the hip. Each hip is scanned. Up to three attempts are permitted to establish correct placement.
- *Knee:* A high resolution image of each knee is performed.
- *Lateral spine:* A high resolution image of the spine (from L4 to T4) in the lateral plane is acquired with the participant lying on their side.
- *Analysis:* Scans at the whole body, lumbar spine and hip sites are analysed by the radiographer at, or soon after, acquisition to generate all numerical measures of bone mass and body composition. Images (whole body, hip, knee, AP lumbar spine and lateral thoraco-lumbar spine) are provided without further analysis for future analysis by researchers.

5. Download/transfer of data

Numerical measures of bone area, bone mineral content and bone mineral density, together with measures of lean and fat mass (for whole body scans only), are directly transferred from the DXA instrument to UK Biobank servers. These data require no post-processing. High resolution images of hips, knees, whole body, AP lumbar spine and lateral thoraco-lumbar spine are exported as Dicom files.

6. QC protocol

All operators undergo standardised central training and scans are performed according to standard operating procedures. Every scan acquisition is assessed in real time by the operator for completeness, movement artefact and presence of any foreign bodies. Findings are recorded and amalgamated into the dataset. Each scan is analysed with optimisation of markers. As part of the pilot phase the first 2000 participant scans were assessed by the acquiring radiographer and by a single experienced senior operator. The intraclass correlation coefficient for measures such as hip bone mineral density was 0.99 demonstrating excellent agreement. A similar “gold standard” assessment of 50 participant scans is undertaken monthly by the same single experienced senior operator and any discrepancies in artefact reporting, analysis or patient positioning are fed back to the technician staff, with retraining if necessary.

7. Data collected

The following data were collected and are available in Showcase:

Bone size, mineral and density by DXA (71 data-fields):

| Field ID | Description |
|----------|---|
| 23220 | Arm BMC (bone mineral content) (left) |
| 23305 | Head bone area |
| 23222 | Arm BMC (bone mineral content) (right) |
| 23317 | Arms combined bone area |
| 23221 | Arm BMD (bone mineral density) (left) |
| 23313 | Arm bone area (left) |
| 23223 | Arm BMD (bone mineral density) (right) |
| 23314 | Arm bone area (right) |
| 23224 | Arms BMC (bone mineral content) |
| 23309 | Ribs bone area |
| 23225 | Arms BMD (bone mineral density) |
| 23311 | Spine bone area |
| 23302 | Femur lower neck BMD (bone mineral density) (left) |
| 23304 | Trunk bone area |
| 23206 | Femur lower neck BMD (bone mineral density) (right) |
| 23307 | Pelvis bone area |

23294 Femur lower neck BMD (bone mineral density) T-score (left)
23318 Legs combined bone area
23207 Femur lower neck BMD (bone mineral density) T-score (right)
23315 Leg bone area (left)
23299 Femur neck BMD (bone mineral density) (left)
23316 Leg bone area (right)
23208 Femur neck BMD (bone mineral density) (right)
23300 Femur neck BMD (bone mineral density) T-score (left)
23209 Femur neck BMD (bone mineral density) T-score (right)
23290 Femur shaft BMD (bone mineral density) (left)
23210 Femur shaft BMD (bone mineral density) (right)
23303 Femur shaft BMD (bone mineral density) T-score (left)
23211 Femur shaft BMD (bone mineral density) T-score (right)
23291 Femur total BMD (bone mineral density) (left)
23212 Femur total BMD (bone mineral density) (right)
23293 Femur total BMD (bone mineral density) T-score (left)
23213 Femur total BMD (bone mineral density) T-score (right)
23295 Femur troch BMD (bone mineral density) (left)
23214 Femur troch BMD (bone mineral density) (right)
23298 Femur troch BMD (bone mineral density) T-score (left)
23215 Femur troch BMD (bone mineral density) T-score (right)
23292 Femur upper neck BMD (bone mineral density) (left)
23216 Femur upper neck BMD (bone mineral density) (right)
23296 Femur upper neck BMD (bone mineral density) T-score (left)
23217 Femur upper neck BMD (bone mineral density) T-score (right)
23297 Femur wards BMD (bone mineral density) (left)
23218 Femur wards BMD (bone mineral density) (right)
23301 Femur wards BMD (bone mineral density) T-score (left)
23219 Femur wards BMD (bone mineral density) T-score (right)
23306 Head BMC (bone mineral content)
23226 Head BMD (bone mineral density)
23203 L1-L4 BMC (bone mineral content)
23204 L1-L4 BMD (bone mineral density)
23205 L1-L4 BMD (bone mineral density) T-score
23320 Leg BMC (bone mineral content) (left)
23228 Leg BMC (bone mineral content) (right)
23227 Leg BMD (bone mineral density) (left)
23229 Leg BMD (bone mineral density) (right)
23230 Legs BMC (bone mineral content)
23231 Legs BMD (bone mineral density)
23308 Pelvis BMC (bone mineral content)
23232 Pelvis BMD (bone mineral density)
23310 Ribs BMC (bone mineral content)
23233 Ribs BMD (bone mineral density)
23312 Spine BMC (bone mineral content)
23234 Spine BMD (bone mineral density)

| | |
|-------|--|
| 23235 | Total BMC (bone mineral content) |
| 23236 | Total BMD (bone mineral density) |
| 23237 | Total BMD (bone mineral density) (left) |
| 23238 | Total BMD (bone mineral density) (right) |
| 23239 | Total BMD (bone mineral density) T-score |
| 23240 | Trunk BMC (bone mineral content) |
| 23241 | Trunk BMD (bone mineral density) |
| 23242 | Trunk BMD (bone mineral density) (left) |
| 23243 | Trunk BMD (bone mineral density) (right) |

Body composition by DXA (49 data-fields):

| Field ID | Description |
|-----------------|-----------------------------------|
| 23244 | Android bone mass |
| 23245 | Android fat mass |
| 23246 | Android lean mass |
| 23247 | Android tissue fat percentage |
| 23248 | Android total mass |
| 23249 | Arm fat mass (left) |
| 23253 | Arm fat mass (right) |
| 23250 | Arm lean mass (left) |
| 23254 | Arm lean mass (right) |
| 23251 | Arm tissue fat percentage (left) |
| 23255 | Arm tissue fat percentage (right) |
| 23252 | Arm total mass (left) |
| 23256 | Arm total mass (right) |
| 23257 | Arms fat mass |
| 23258 | Arms lean mass |
| 23259 | Arms tissue fat percentage |
| 23260 | Arms total mass |
| 23261 | Gynoid bone mass |
| 23262 | Gynoid fat mass |
| 23263 | Gynoid lean mass |
| 23264 | Gynoid tissue fat percentage |
| 23265 | Gynoid total mass |
| 23200 | L1-L4 area |
| 23201 | L1-L4 average height |
| 23202 | L1-L4 average width |
| 23266 | Leg fat mass (left) |
| 23270 | Leg fat mass (right) |
| 23267 | Leg lean mass (left) |
| 23271 | Leg lean mass (right) |
| 23268 | Leg tissue fat percentage (left) |
| 23272 | Leg tissue fat percentage (right) |
| 23269 | Leg total mass (left) |
| 23273 | Leg total mass (right) |

| | |
|-------|--------------------------------------|
| 23274 | Legs fat mass |
| 23275 | Legs lean mass |
| 23276 | Legs tissue fat percentage |
| 23277 | Legs total mass |
| 23278 | Total fat mass |
| 23279 | Total fat-free mass |
| 23280 | Total lean mass |
| 23281 | Total tissue fat percentage |
| 23282 | Total tissue mass |
| 23283 | Total total mass |
| 23284 | Trunk fat mass |
| 23285 | Trunk lean mass |
| 23286 | Trunk tissue fat percentage |
| 23287 | Trunk total mass |
| 23288 | VAT (visceral adipose tissue) mass |
| 23289 | VAT (visceral adipose tissue) volume |
